
Combure® Test

REF 11896962 257

cobas
▼ 50

English

Intended use

Six-patch test strip for the semi-quantitative determination of glucose, leukocytes, nitrite, protein, urobilinogen and blood in urine by visual reading.

For professional use only.

Summary

Urine test strips are used to measure certain constituents in urine which are significant of renal, urinary, hepatic and metabolic disorders.

Test principle

Leukocytes (LEU): The test reveals the presence of granulocyte esterases. These esterases cleave an indoxyl ester, and the indoxyl so liberated reacts with a diazonium salt to produce a violet dye. Bacteria, trichomonads or erythrocytes present in the urine do not affect the reaction.

Nitrite (NIT): The test is based on the principle of the Griess test and is specific for nitrite. The reaction reveals the presence of nitrite and hence indirectly nitrite-forming bacteria in the urine by a pink-to-red coloration of the test patch. Even a slight pink coloration is indicative of significant bacteruria.

Protein (PRO): The test is based on the principle of the protein error of a pH indicator. It is particularly sensitive to albumin. Quinine, quinidine, chloroquine, tolbutamide and an elevated pH (up to 9) do not affect the test.

Glucose (GLU): The glucose determination is based on the specific glucose-oxidase/peroxidase reaction (GOD/POD method). The test is independent of the pH and specific gravity of the urine and is not affected by the presence of ketone bodies.

Urobilinogen (UBG): A stable diazonium salt reacts almost immediately with urobilinogen to give a red azo dye. The test is specific for urobilinogen and is not susceptible to the interfering factors known to affect the Ehrlich's test.

Blood (ERY/Hb): The peroxidase-like action of hemoglobin and myoglobin specifically catalyzes the oxidation of the indicator by means of the organic hydroperoxide contained in the test paper to give a blue-green coloration.

Reagents

Each test contains per 1 cm² test patch area the following:

Leukocytes: Indoxylcarboxylic acid ester 15.5 µg; methoxymorpholinobenzene diazonium salt 5.5 µg

Nitrite: 3-hydroxy-1,2,3,4-tetrahydro-7,8-benzoquinoline 33.5 µg; sulfanilamide 29.1 µg

Protein: 3',3",5'-tetrachlorophenol-3,4,5,6-tetrabromosulfophthalein 13.9 µg

Glucose: 3',3",5'-tetramethylbenzidine 103.5 µg; GOD 6 U, POD 35 U

Urobilinogen: 4-methoxybenzene-diazonium-tetrafluoroborate 67.7 µg

Blood: 3',3",5'-tetramethylbenzidine 52.8 µg; 2,5-dimethyl-2,5-dihydroperoxyhexane 297.2 µg

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines.

Safety data sheet available for professional user on request.

The stopper of the test strip vial contains a non-toxic silicate-based desiccant which must not be removed. If ingested by accident, drink large quantities of water.

Reagent handling

Test strips are ready for use.

Storage and stability

Store the package at 2-30 °C. The test strips are stable up to the expiration date specified on the box, when stored in the original container.

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Use fresh urine that has not been centrifuged. The urine specimen should not stand for more than 2 hours before testing. In case of longer standing, mix before use.

Use only clean, well-rinsed vessels to collect urine.

Do not add preservatives to the urine.

Materials provided

REF 118969622257, package with 50 test strips

Materials required (but not provided)

General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document.

1. Use fresh urine that has not been centrifuged. Thoroughly mix the urine sample. The sample should be at room temperature when the test is performed and should not have been standing for more than 2 hours.

2. Take a test strip out of the container. Close the container again with the original desiccant stopper immediately after removal of the strip. This is important as otherwise the test areas may become discolored due to environmental influences such as moisture or nitrite gases in the air. Incorrect results may be obtained.

3. Briefly (about 1 second) dip the test strip into the urine making sure that all test areas are moistened.

4. When withdrawing the test strip, wipe the edge against the rim of the vessel to remove excess urine.

5. After 60 seconds (60-120 seconds for the leukocyte test area), compare the reaction colors of the test areas with the colors on the label and assign always the value of the nearest color block.

- Compare the 6th (blood) test area with both color scales as separate color scales are given for erythrocytes and hemoglobin.

Any color changes appearing only along the edges of the test areas, or developing after more than 2 minutes, do not have any diagnostic significance.

Quality control

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

For quality control, use commercially available urine controls, or other suitable control material. Important note for reporting results (for professional users).

According to the regulations from the German Medical Association for quality assurance of medical laboratory analyses dated 11/23/2007, the decision to classify a laboratory test result to either part B1 or B2 depends on the way the test results are expressed in the report (scale level).

The specification in the report defines whether a determination is quantitative or qualitative and therefore which legal requirements for quality assurance (B1 for quantitative or B2 for qualitative) need to be followed. Examples of qualitative characteristics include titer levels, concentrations/color ranges (+ to ++++) or a defined range of values. A characteristic of a quantitative value is when the value has a corresponding measured unit value.

Limitations - interference

Leukocytes: Formaldeyde (stabilizer) and medication with imipenem, meropenem and clavulanic acid may cause false-positive reactions. If the urine specimen has a pronounced intrinsic color (for example due to the presence of bilirubin or nitrofurantoin), the reaction color may be intensified due to an additive effect. Urinary protein excretions in excess of 500 mg/dL and urinary glucose excretions in excess of 1 g/dL may diminish the intensity of the reaction color, as can cephalaxin and gentamicin if administered in high daily doses, or boric acid if used as a preservative.

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- Accuracy**
- Leukozyten: ≥ 90 % compared with counting chamber, Nitrite: ≥ 90 % for 10⁷ gram-negative organisms compared with Griess test, Protein: 90 % compared with radial immuno-diffusion, Glucose: ≥ 90 % compared with hexokinase method, Urobilinogen: ≥ 95 % compared with Watson & Henry method, Blood: ≥ 90 % compared with counting chamber.

- Estimated values**
- Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

- For further information, please refer to the appropriate operator's manual for the instrument concerned, and the Method Sheets of all necessary components.
- A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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