
Combur⁹ Test

[REF] 04510038 191

[REF] 04510046 040

English
Intended use
Nine-patch test strip for the semi-quantitative determination of leukocytes, nitrite, pH, protein, glucose, ketone bodies, urobilinogen, bilirubin and blood in urine by visual reading.
For professional use only.

Summary
Urine test strips are used to measure certain constituents in urine which are significant of renal, urinary, hepatic and metabolic disorders.

Test principle

Leukocytes (LEU): The test reveals the presence of granulocyte esterases. These esterases cleave an indoxyl ester, and the indoxyl so liberated reacts with a diazonium salt to produce a violet dye. Bacteria, trichomonads or erythrocytes present in the urine do not affect the reaction.

Nitrite (NIT): The test is based on the principle of the Griess test and is specific for nitrite. The reaction reveals the presence of nitrite and hence indirectly nitrite-forming bacteria in the urine by a pink-to-red coloration of the test patch. Even a slight pink coloration is indicative of significant bacteriuria.

pH: The test paper contains the indicators methyl red, phenolphthalein and bromothymol blue and reacts specifically with H⁺-ions. The most frequent pH values of fresh urine from healthy subjects lie between 5 and 6.

Protein (PRO): The test is based on the principle of the protein error of a pH indicator. It is particularly sensitive to albumin. An elevated pH (up to 9) does not affect the test.

Glucose (GLU): The glucose determination is based on the specific glucose-oxidase/oxidase reaction (GOD/POD method). The test is independent of the pH and specific gravity of the urine and is not affected by the presence of ketone bodies.

Ketones (KET): This test is based on the principle of Legal’s test and is more sensitive to acetoacetic acid than to acetone.

Urobilinogen (UBG): A stable diazonium salt reacts almost immediately with urobilinogen to give a red azo dye. The test is specific for urobilinogen and is not susceptible to the interfering factors known to affect the Ehrlich’s test.

Bilirubin (BL): The test is based on the coupling of bilirubin with a diazonium salt. Even the slightest pink coloration constitutes a positive, i.e. pathologic, result. Other urinary constituents produce a more or less intense yellow coloration.

Blood (ERY/Hb): The peroxidase-like action of hemoglobin and myoglobin specifically catalyzes the oxidation of the indicator by means of the organic hydroperoxide contained in the test paper to give a blue-green coloration.

Reagents

Each test contains per 1 cm² test patch area the following:

Leukocytes: Indoxylcarbonic acid ester 15.5 µg; methoxymorpholinobenzene diazonium salt 5.5 µg

Nitrite: 3-hydroxy-1,2,3,4-tetrahydro-7,8-benzoquinoline 33.5 µg; sulfanilamide 29.1 µg

pH: Bromothymol blue 13.9 µg; methyl red 1.2 µg; phenolphthalein 8.6 µg

Protein: 3',3',5',5"-tetrachlorophenol-3,4,5,6-tetrabromsulphthalein 13.9 µg

Glucose: 3,3',5',5"-tetramethylbenzidine 103.5 µg; GOD 6 U, POD 35 U

Ketones: Sodium nitroprusside 157.2 µg

Urobilinogen: 4-methoxybenzene-diazonium-tetrafluoroborate 67.7 µg

Bilirubin: 2,6-dichlorobenzene-diazonium-tetrafluoroborate 16.7 µg

Blood: 3,3',5',5"-tetramethylbenzidine 52.8 µg; 2,5-dimethyl-2,5-dihydroperoxyhexane 297.2 µg

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

The stopper of the test strip vial contains a non-toxic silicate-based desiccant which must not be removed. If ingested by accident, drink large quantities of water.

Reagent handling

Test strips are ready for use.

Storage and stability

Store the package at 2-30 °C. The test strips are stable up to the expiration date specified on the box, when stored in the original container.

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

Specimen collection and preparation

Use only a fresh midstream urine sample¹ that has not been centrifuged.

For specimen collection and preparation only use suitable tubes or collection containers.

Use only clean, well-rinsed vessels to collect urine.

False-positive readings for blood and glucose can result from residues of strongly oxidizing disinfectants in the specimen collection vessel. The urine specimen should not stand for more than 2 hours before testing. In case of longer standing, store the specimen refrigerated at 4° C. At the time of testing, ensure the sample is at room temperature and mix the sample thoroughly before use. Longer standing times can lead to false results owing to the following influences: Proliferation of bacteria; A rise in pH due to ammonia formed as a result of bacterial degradation of urea.

Do not add preservatives to the urine.

Do not expose urine specimens to sunlight as this induces oxidation of bilirubin and urobilinogen and hence leads to artificially low results for these two parameters.² Drugs that turn red in an acid environment (e.g. phenazopyridin) may produce false- positive readings or reddish colorations on the test patches for nitrite, protein, urobilinogen and bilirubin.

Materials provided

- [REF] 04510038191, package with 50 test strips

- [REF] 04510046040, package with 100 test strips

Materials required (but not provided)

- General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document.

- Use fresh urine that has not been centrifuged. Thoroughly mix the urine sample. The sample should be at room temperature when the test is performed and should not have been standing for more than 2 hours.
- Take a test strip out of the container. Close the container again with the original desiccant stopper immediately after removal of the strip. This is important as otherwise the test areas may become discolored due to environmental influences such as moisture or nitrite gases in the air. Incorrect results may be obtained.
- Briefly (about 1 second) dip the test strip into the urine making sure that all test areas are moistened.
- When withdrawing the test strip, wipe the edge against the rim of the vessel to remove excess urine.
- After 60 seconds (60-120 seconds for the leukocyte test area), compare the reaction colors of the test areas with the colors on the label and assign always the value of the nearest color block.

- Compare the 9th (blood) test area with both color scales as separate color scales are given for erythrocytes and hemoglobin.

Any color changes appearing only along the edges of the test areas, or developing after more than 2 minutes, do not have any diagnostic significance.

Quality control

The control intervals and limits should be adapted to each laboratory’s individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

For quality control, use commercially available urine controls, or other suitable control material.

Important note for reporting results (for professional users). According to the regulations from the German Medical Association for quality assurance of medical laboratory analyses dated 11/23/2007, the decision to classify a laboratory test result to either part B1 or B2 depends on the way the test results are expressed in the report (scale level).

The specification in the report defines whether a determination is quantitative or qualitative and therefore which legal requirements for quality assurance (B1 for quantitative or B2 for qualitative) need to be followed. Examples of qualitative characteristics include liter levels, concentrations/color ranges (+ to +++) or a defined range of values. A characteristic of a quantitative value is when the value has a corresponding measured unit value.

Limitations - interference

Leukocytes: Formaldehyde (stabilizer) and medication with imipenem, meropenem and clavulanic acid may cause false-positive reactions. If the urine specimen has a pronounced intrinsic color (for example due to the presence of bilirubin or nitrofurantoin), the reaction color may be intensified due to an additive effect. Urinary protein excretions in excess of 500 mg/dL and urinary glucose excretions in excess of 1 g/dL may diminish the intensity of the reaction color, as can cephalaxin and gentamicin if administered in high daily doses, or boric acid if used as a preservative.

Vaginal secretion or menstrual blood may contaminate urine from females. This may be misleading and can be minimized by tamponing the vagina when acute symptoms necessitate examination of urine¹.

Nitrite:Prolonged urinary retention in the bladder (4-8 hours) is essential in order to obtain an accurate result. Administration of antibiotics or chemical drugs should be discontinued 3 days before the test. 95 % of all bacteria responsible for urinary tract infections are Gram-negative rods (E.coli, Proteus mirabilis, Klebsiella pneumoniae), present singular or in short chains, usually in large numbers. In hospitalized patients Gram-positive cocci (enterococci) are of significance². Normal nutrition as a rule ensures a sufficiently high content of nitrate in the urine for the detection of bacteria. The most frequent pathogen responsible for urinary tract infections, E. coli and most of the other urine-pathogenic organisms (Klebsiella, Aerobacter, Citrobacter, Salmonella and to some extent also enterococci, staphylococci, and Pseudomonas) reduce urinary nitrate to nitrite and can therefore be detected indirectly with test strips³. Some common uropathogens, e.g. Enterococcus spp. and Staphylococcus spp., do not reduce urinary nitrate to nitrite and will therefore not be detected whatever their urinary concentration¹. False-negative results may occur as a result of strong diuresis with frequent voiding of urine³, insufficient intake or too short retention of urine in the bladder². Large amounts of ascorbic acid decrease the sensitivity of the test.*Attention:* Nitrogen oxides present in the atmosphere may have an influence on the stability of the nitrite test pad¹.

Protein: False-positive readings may be found after infusion of polyvinylpyrrolidone (blood substitute), or when the urine collection vessel contains chlorhexidine or traces of disinfectants possessing quaternary ammonium groups.

Glucose: The effect of ascorbic acid has been largely eliminated so that at glucose concentrations of 100 mg/dL and above even high ascorbic acid concentrations are not likely to give false-negative results.

Ketones: Phenylketones and phthalein compounds produce red colors on the test patch; they are, however, quite distinct from the violet colors produced by ketone bodies and can lead to false-positive results. Captopril, mesna (2-Mercaptoethansulfonic acid sodium salt) and other substances containing sulhydryl groups may produce false-positive results.

Urobilinogen: Nitrite concentrations above 5 mg/dL or formaldehyde (stabilizer) above 200 mg/dL may cause a decrease in the color reaction.

Bilirubin: Large amounts of ascorbic acid lower the sensitivity of the test.

Do not expose urine specimens to sunlight as this induces oxidation of bilirubin and urobilinogen and hence leads to artificially low results for these two parameters.

Drugs that turn red in an acidic environment (e.g. phenazopyridine) may produce false-positive readings or reddish colorations on the test patches for nitrite, protein, urobilinogen and bilirubin.

Blood: The values appearing on the printout refer to intact erythrocytes. At concentrations of about 5-50 Ery/µL, significant hemolysis (such as may occur on prolonged standing of the urine) leads to values which are higher than the corresponding concentrations given for intact erythrocytes. Ascorbic acid has virtually no effect on the test. In women the test for blood may be falsified from 3 days before to 3 days after a period. It is therefore advisable not to perform the test during this time. After physical activity, e.g. strenuous jogging, raised values for erythrocytes and protein may occur without being signs of disease.

False-positive readings, particularly for erythrocytes, glucose and protein can result from residues of strongly oxidizing disinfectants in the specimen collection vessel.

Note:

Knowledge of the effects of drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is therefore advisable to repeat the test after discontinuing a particular drug.

For diagnostic purposes, the results should always be assessed in conjunction with the patient’s medical history, clinical examination and other findings.

Limits and ranges

By visual reading

• **Measuring range**

Leukocytes: NEG – – 500 LEU/µL (3+), **Nitrite:** NEG – POS, **pH:** 5-9, **Protein:** NEG - 500 mg/dL (5 g/L; 3+), **Glucose:** NORM – 1000 mg/dL (55 mmol/L; 4+), **Ketone bodies:** NEG – 150 mg/dL (15 mmol/L; 3+), **Urobilinogen:** NORM – 12 mg/dL (200 µmol/L; 4+), **Bilirubin:** NEG – – 6 mg/dL (100 µmol/L; 3+), **Blood:** NEG – – 250 ERY/µL (4+).

• **Lower limits of measurement**

Lower detection limit

Leukocytes: 10-25 LEU/µL, **Nitrite:** 0.05 mg/dL (11 µmol/L), **Protein:** 6 mg albumin/dL, **Glucose:** 40 mg/dL (2.2 mmol/L), **Ketone bodies:** For aceto-acetic acid 5 mg/dL (0.5 mmol/L), **Urobilinogen:** 0.4 mg/dL (7 µmol/L), **Bilirubin:** 0.5 mg/dL (9 µmol/L), **Blood:** intact erythrocytes: 5 ERY/µL, Hemoglobin or hemolysed erythrocytes: corresponding to 10 ERY/µL.

• **Accuracy**

Leukocytes: ≥ 90 % compared with counting chamber, **Nitrite:** ≥ 90 % for 10⁷ gram-negative organisms compared with Griess test, **pH:** ≥ 95 % compared with pH-meter, **Protein:** 90 % compared with radial immuno-diffusion, **Glucose:** ≥ 90 % compared with hexokinase method, **Ketone bodies:** ≥ 85 % compared with photometric enzymatic determination of acetate, **Urobilinogen:** ≥ 95 % compared with Watson & Henry method, **Bilirubin:** ≥ 85 % compared with total bilirubin determination by Jendrassik’s method (direct bilirubin), **Blood:** ≥ 90 % compared with counting chamber.

Expected values

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

For further information, please refer to the appropriate operator’s manual for the instrument concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

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- ↑ T. Kouri, G. Fogazzi, V. Gant et al. European Urinalysis Guidelines. Scandinavian Journal of Clinical & Laboratory Investigation 2000; Vol 60: Suppl 231.
- ↑ REF 12254320001 Compendium of Urinalysis, Urine test strips and microscopy
- ↑ REF 12254620001 Compendium of Urinalysis, Urinalysis with test strips

Symbols / Symbole / Símbolos

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard. / In Erweiterung zur ISO 15223-1 werden von Roche Diagnostics folgende Symbole und Zeichen verwendet. / Roche Diagnostics emplea los siguientes símbolos y signos adicionalmente a los indicados en la norma ISO 15223-1.

	Global Trade Item Number / Globale Artikelnummer GTIN / Número mundial de artículo comercial
	

COBAS and COMBUR-TEST are trademarks of Roche.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2015, Roche Diagnostics

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim

www.roche.com