

Influenza A/B+COVID-19/RSV Combo Ag Test (self-testing)

Influenza A/B+COVID-19/RSV Combo Ag Test is an in vitro immunochromatographic assay for the qualitative and differential detection of nucleocapsid protein antigen from influenza A (including the subtype H1N1), influenza B, respiratory syncytial virus and/or SARS-CoV-2 in nasal swab specimens. It is intended to aid in the rapid diagnosis of influenza A, influenza B, respiratory syncytial virus and/or SARS-CoV-2 infections. This test provides only a preliminary test result. Therefore, any reactive specimen with the Influenza A/B+COVID-19/RSV Combo Ag Test must be confirmed with alternative testing method(s) and clinical findings.

Advantage:

- ★ Non-invasive sampling
- ★ Rapid Test (Test result available in 15 min.)
- ★ Read test results visually.
- ★ Do not require test equipment.
- ★ Simple Operation.




Certification





ISO13485



Core Technology Co.,Ltd.

Chinese Address: No.30 Area 9,Douda Avenue,
Fangshan District,Beijing,102433,China
Tel:86-10-69390616 Website:www.coretests.com
Email: info@coretests.com

CorDx Inc.


USA Address:9540 Waples St,Unit C, San Diego,
CA 92121 USA
Tel:+1 (858) 3331122 Website:www.coretests.us
Email: info@coretests.com



Scan or Click QR Code

Influenza A/B+COVID-19/RSV Combo Ag Test (self-testing)

Test Procedure (Nasal Swab Specimen)



1 remove foil
Remove the foil from the top of the sample collection tube.

2
Place the tube in the tube rack.

3 not touch swab head
Remove a sampling swab from the pouch.

4 insert 2-4 cm
Using the sterile swab provided in the kit, carefully insert the swab into one nostril.

5 roll 5 times
roll 5 times
The swab tip should be inserted up to 2-4 cm until resistance is met. Roll the swab 5 times in a circular motion around the inside wall to ensure that both mucus and cells are collected. Using the same swab, repeat this process for the other nostril to ensure that an adequate sample is collected from both nasal cavities.

6
Withdraw the swab from the nasal cavity.

7 rotate 10 times
The specimen is now ready for preparation using the extraction buffer provided in the test kit. Insert the swab in collection tube to the bottom, rotate and squeeze the swab 10 times while pressing the head against the bottom and side of the collection tube.

8 wait 1 min
Leave the swab in the collection tube for 1 minute.

9 squeeze
Rotate and squeeze the tube several times with fingers from outside of the tube to immerse the swab. Remove the swab.

10
Attach the dropper tip firmly onto the tube. Mix thoroughly by swirling or flicking the bottom of the tube.

11
Add 3 drops of specimen to the sample well.

12
Read results within 15 minutes.

13 Interpretation of Results

For Flu A/B


(+)	(+)	(+)	(-)
C T2 T1	C T2 T1	C T2 T1	C T2 T1
Flu A/B	Flu B	Flu A	Neg

Invalid

For COVID-19/RSV

(+)	(+)	(+)	(-)
C T2 T1	C T2 T1	C T2 T1	C T2 T1
COVID-19	COVID-19	RSV	Neg

Invalid

Product	Brand	Specimen	Format	REF Numbers	Package Size	Certificate
Influenza A/B+ COVID-19/RSV Combo Ag Test	CorDx	Nasal Swab	Cassette	BC292-01	1 T	
				BC292-02	2 Ts	
				BC292-05	5 Ts	
				BC292-25	25 Ts	

Core Technology Co.,Ltd.

Chinese Address: No.30 Area 9,Douda Avenue,
Fangshan District,Beijing,102433,China
Tel:86-10-69390616 Website:www.coretests.com
Email: info@coretests.com

CorDx Inc.

USA Address:9540 Waples St.Unit C, San Diego,
CA 92121 USA
Tel:+1 (858)3331122 Website:www.coretests.us
Email: info@coretests.com



Scan or Click QR Code