

Technical Data (TD)

Urine Container 100 ml w/Screw Cap				
REF	06-065-1004	06-065-1104		

1.	Description / Specification	Cup	Cap	
1.01	Description	Universal Container 100 ml w/Screw Cap		
1.02	Dimension	See Customer Drawing		
1.03	Volume	100 ml		
1.04	Material / Resin	Polypropylene (PP)	Polyethylene (HDPE)	
	Grade of Resin	MAT123	MAT131	
	Properties of Resin	See www.boettger.com		
1.05	REF	06-065-1*04		
	Colour Cap	Green		
	Grade of Pigment	PIG114		
1.06	Hygiene Level	Standard Hygienic		
1.07	Quality Control	Testing of attributive and variable characteristics in accordance with the valid specification		
1.08	Storage	To be stored in closed original packaging: Dry, at room temperature		
1.09	Usage / Purpose	General laboratory product, Collection and evaluation of samples, also for urine		
1.10	Conformities / Regulations	See Classification		
1.11	Product Stewardship	See Product Stewardship of Grades of Resin/Pigments		
1.12	Other information	For single use only		
2.	Features	Cup	Cap	
2.01	Basic features	-		
2.02	Temperature range	-196°C to 121°C	-50°C to 70°C	
2.03	Autoclavability	Yes	No	
2.04	Centrifugation	N/A		
2.05	Chemical Resistance	See www.boettger.com		
2.06	Expiry	5 years from LOT (YYYY-MM)		
2.07	Other information	-		
3.	Packaging	06-065-1004	06-065-1104	
3.01	Container Pcs/Bag	25	25	
	Caps Pcs/Bag	250	100	
3.02	Pcs/Ctn	500	500	
3.03	Pcs/Pallet	12000	12000	
3.06	LOT / VP	WWYY (week, year of production) / XX (packaging)		
3.07	Other information	See Measures		
4.	Other Information			
4.01	Codes and Numbers	See Codes and Numbers		
4.02	Symbols	See Symbols		






All informations without guarantee. Above mentioned products are manufactured according to company specifications and our certified Quality Management System corresponding to DIN EN ISO 13485. The quality check is carried out on a sample basis; an individual examination is not possible. Deficiencies parts may occur. We accept no liability for consequential damages.

Revision	Drawn	Approved	Confidential: Information contained in this document or drawing is confidential and proprietary to Paul Boettger GmbH & Co. KG. This document may not be reproduced for any reason without written permission from Paul Boettger GmbH & Co. KG. All rights of design, invention and copyright are reserved.
10	22.07.2022 QMB	22.07.2022 QMB	



boettger

Technical Data (TD)

5.	Label (Example) according to IVDR 2017/746/EC
<div style="display: flex; justify-content: space-between;"> <div style="width: 60%;"> <p>Paul Boettger GmbH & Co. KG Miesleuthenweg 35 - 45 94249 Bodenmais/Germany www.boettger.com</p> </div> <div style="width: 35%; text-align: right;">  boettger </div> </div> <hr/> <div style="display: flex; justify-content: space-between; align-items: center;"> <div style="font-size: 24px; font-weight: bold;">500</div> <div style="font-size: 24px;">St/pcs</div> <div style="border: 1px solid black; padding: 2px 5px; font-weight: bold;">REF</div> <div style="font-size: 24px; font-weight: bold;">06-065-1004</div> </div> <hr/> <p style="font-size: 18px;"> Urinbecher SV 100 m/Deckel (grün) Urine Cup 100 PP w/screw cap (green) </p> <div style="text-align: right; font-size: 10px; margin-top: 10px;">V210601</div> <hr/> <div style="display: flex; justify-content: space-between; align-items: center;"> <div style="width: 60%;">  <div style="display: flex; align-items: center;"> <div style="border: 1px solid black; padding: 2px 5px; font-weight: bold; margin-right: 5px;">IVD</div> <div style="font-size: 10px; margin-right: 5px;"> Paul Boettger GmbH Co. KG D-94249 Bodenmais </div> <div style="font-size: 10px; margin-right: 5px;">  2026-06 </div> </div> </div> <div style="width: 35%; text-align: right;">  </div> </div> <p style="font-size: 10px; margin-top: 5px;">Lagern - trocken bei Raumtemperatur / Store - dry at room temperature</p> <hr/> <div style="display: flex; justify-content: space-between; align-items: center;"> <div style="width: 30%;"> <div style="border: 1px solid black; padding: 2px 5px; font-weight: bold; margin-bottom: 5px;">UDI</div>  <p style="font-size: 8px;">(01) 0 4042836 00184 0 (17) 260600 (10) 2321</p> </div> <div style="width: 65%; text-align: right;"> <div style="font-size: 24px; font-weight: bold; margin-bottom: 5px;">2321</div> <div style="border: 1px solid black; padding: 2px 5px; font-weight: bold; margin-right: 5px;">LOT</div> <div style="font-size: 24px; font-weight: bold; margin-bottom: 5px;">12</div> <div style="font-size: 24px; font-weight: bold; margin-bottom: 5px;">VP</div> <div style="font-size: 18px;">Ctn 14</div> </div> </div>	

Revision	Drawn	Approved	Confidential: Information contained in this document or drawing is confidential and proprietary to Paul Boettger GmbH & Co. KG. This document may not be reproduced for any reason without written permission from Paul Boettger GmbH & Co. KG. All rights of design, invention and copyright are reserved.
10	22.07.2022 QMB	22.07.2022 QMB	